

Xyltech™ Cell Proliferation Control System

Example




Fig.1 Changes in the number of normal hADSCs.

hADSCs were cultured using Growth MSC+Supplement XF [for proliferation] (Growth MSC(XF)) and MSC-01 Xeno-Free [for proliferation control] (MSC-01). When the medium was replaced on day 3, it was confirmed that cells quickly repopulated upon switching from MSC-01 to Growth MSC (XF).




Fig.2 Phase contrast microscope images of normal hADSCs.

Features

- Xyltech™ Growth MSC (with Supplement XF or Supplement AF) enables the culture of human adipose-derived mesenchymal stem cells without the need for a coating agent.
- Xyltech™ MSC-01 Xeno-Free and Xyltech™ MSC-02 Animal-Free can be cultured maintaining cell properties while suppressing proliferation.
- Xyltech™ MSC series proposes various culture applications.

Cell proliferation control

Reducing the number of passages and managing experimental schedules

Cell freezing avoidance

Short-term cell storage and avoidance of cultivation tasks during holidays

Non-freezing transport of cells

Cell density control during live cell transport

ブルボン再生医科学研究所

Bourbon Biomedical Advanced Research Laboratories, Inc.
1-3-1, Ekimae, Kashiwazaki City, Niigata Pref.,
945-0055 JAPAN
E-mail support@bourbon-barl.co.jp
URL <https://www.bourbon-barl.co.jp/eg/>



Serum-Free hMSC(human mesenchymal stem cell)Proliferation Control Medium

Xyltech™ MSC-01 Xeno-Free

Xyltech™ MSC-02 Animal-Free

Serum-Free hMSC Growth Medium

Xyltech™ Growth MSC

Products

Cat. No.	Product Name	Features / Applications	Volume
10401	Xyltech™ MSC-01 Xeno-Free	Serum-free culture medium for human Mesenchymal Stem Cells [Cell proliferation suppression] (Xeno-Free)	100ml(P)
10402	Xyltech™ MSC-02 Animal-Free	Serum-free culture medium for human Mesenchymal Stem Cells [Cell proliferation suppression] (Animal-Free)	100ml(P)
10411	Xyltech™ Growth MSC Medium	Serum-free culture medium for human Mesenchymal Stem Cells [Cell proliferation] Basal culture medium for MSCs	500ml(P)
10412	Xyltech™ Growth MSC Supplement XF	Serum-free culture medium for human Mesenchymal Stem Cells [Cell proliferation] Supplement (Xeno-Free)	10ml(P)
10413	Xyltech™ Growth MSC Supplement AF	Serum-free culture medium for human Mesenchymal Stem Cells [Cell proliferation] Supplement (Animal-Free)	10ml(P)

(P)PET bottle

These products are research reagents and not intended for human or animal treatment or diagnostic purposes.

Serum-Free hMSC(human mesenchymal stem cell) Proliferation Control Medium

Serum-Free hMSC Growth Medium




Xyltech™ series


Xyltech™ MSC-01 Xeno-Free Xyltech™ MSC-02 Animal-Free

- Human MSCs culture medium for suppressing cell proliferation rate while maintaining the cell properties for up to 3 days
- Rapid recovery of cell proliferation after switching to the growth medium
- Xyltech™ MSC-01 Xeno-Free contains only human-derived components, no animal derivatives. Xyltech™ MSC-02 Animal-Free contains no ingredients derived from animals, including humans.
- These products do not contain glucose

Cell Morphology




Cell Proliferation Profiling



The cell number of normal hADSCs cultured using Growth MSC (for cell proliferation, blue) and MSC-02 Animal-Free (for cell proliferation suppression, orange) was determined.

Cell Differentiation



Differentiation induction of normal hADSCs into adipocytes (left), osteoblasts (middle), and chondrocytes (right). Prior to the differentiation induction, hADSCs were cultured with MSC-01 or MSC-02 (for growth suppression) for 3 days.

We recommend using products Xyltech™ and Xyltech™ Growth as a combination.


Xyltech™ Growth series



Xyltech™ Growth MSC Medium Xyltech™ Growth MSC Supplement XF Xyltech™ Growth MSC Supplement AF


- Human MSCs culture medium for cell proliferation while maintaining the cell properties
- Xyltech™ Growth MSC Medium used as the basal culture medium, supplemented with Xyltech™ Growth MSC Supplement XF or Xyltech™ Growth MSC Supplement AF
- Xyltech™ Growth MSC Supplement XF is xeno-free, Xyltech™ Growth MSC Supplement AF is animal-free
- No coating required on culture vessels

Cell Morphology




Phase contrast microscope images of normal hADSCs after 3 days of culture with Growth MSC.

Cell Differentiation



Differentiation induction of normal hADSCs into adipocytes (left), osteoblasts (middle), and chondrocytes (right). Prior to the differentiation induction, hADSCs were cultured with Growth MSC+Supplement AF medium for 3 days.

Cell Proliferation Profiling



— Growth MSC + Supplement AF — Growth MSC + Supplement XF — Company X (Xeno-Free)

The cell number of normal hADSCs cultured up to the 6th passage using Growth MSC and Company X xeno-free medium was determined.